


SEMINOR RIGID/FLEXIBLE

MULTI-PURPOSE CANNULA FOR ASSISTED HUMAN REPRODUCTION

- A) The Seminor Rigid is a semi-rigid cannula for intrauterine and cervical insemination.
- B) The Seminor Flexible is a cannula for intrauterine and cervical insemination.


FEATURES

- the Seminor Rigid (A) can be shaped according to the degree of uterine flexion
- the Seminor Flexible (B) enables easy insertion with minimal patient discomfort
- the conical proximal end closes the cervix adequately for maximum sperm transfer
- as result of a direct connection between Seminor and syringe, the loss of spermatozoa during the transfer is minimized




GYNO INFO SEMINOR RIGID/FLEXIBLE

Length	48 mm (\pm 0.5 mm)
Volume	min. 0.098 ml
Sterile	sterilized by gamma irradiation
Packaging	125/shipping box
Quality control	MEA and LAL
Order number	A) #4502 A  B) #4502 B 



GMP Belgium Seminor 09/2014



Gynetics Medical Products n.v.
Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
tel. +32 11 64 58 72 / fax +32 11 64 06 77
info@gynetics.be / www.gynetics.be

GYNETICS[®]
Practise better Medicine