

Hyaluronidase

- For removing cumulus cells for oocyte denudation
- Ready to use as adjusted to 80 (IU/mL)
- It can be used up to 7 days after opening



Type	Order Number	Code	Contents
Hyaluronidase	95940	HYBD - 1x5	5x1ml
	95931	HYBD - 10	1x10mL

COMPONENTS

- Calcium Chloride
- Glucose
- HEPES
- Hyaluronidase from Bouine testes
- Magnesium Sulfate
- Potassium Chloride
- Sodium Phosphate
- Human Serum Albumin (4.0g/l)
- Sodium Bicarbonate
- Sodium Chloride
- Lactate
- Pyruvate

QUALITY CONTROL

pH 7.3-7.6 / Osmolality 270-290mOsm/L
 Endotoxin <0.5EU/mL / Mouse Embryo Assay ≥ 80%
 Sterility Test
 Storage: 2-8 °C
 Shelf Life: 12 months

KITAZATO CORPORATION
 100-10 Yanagishima, Fuji
 Shizuoka 416-0932 Japan
 Tel +81 545-65-7122
 Fax +81 545-65-7128
 contact@kitazato.co.jp

KITAZATO | DIBIMED
 C/ Jorge Comin, 3 Bajo 1-2
 46015 Valencia, Spain
 Tel +34 96 305 63 95
 Fax +34 96 305 63 96
 info@dibimed.com
Official Representative of Kitazato

KITAZATO®