

NEW IVF MEDIA

PVP Polyvinylpyrrolidone

- For immobilization of sperm
- Ready to Use
- It can be used up to 7 days after opening



Type	Order Number	Code	Contents
PVP 10%	95910	PVPD10 - 0,2x5	5x0,2ml
	95901	PVPD10 - 1	1x1mL

COMPONENTS

- Calcium Chloride
- Glucose
- HEPES
- Magnesium Sulfate
- 10% Polyvinylpyrrolidone
- Potassium Chloride
- Human Serum Albumin (4.0g/l)
- Sodium Phosphate
- Sodium Bicarbonate
- Sodium Chloride
- Lactate
- Pyruvate

QUALITY CONTROL

pH 7.2-7.6 / Endotoxin <0.5EU/mL / Sterility Test
Osmolality 280-340mOsm/kg
Mouse Embryo Assay ≥ 80%
Human Sperm Survival Assay ≥ 80%
Storage: 2-8 °C
Shelf Life: 9 months

KITA ZATO®